

# Monobloc Fill/Finish Packaging System



MNB-2000  
Monobloc



- Ability to Sort, Feed, Fill, Plug, Stopper, Crimp, Cap, Label, and/or Accumulate in a Compact Footprint
- User Friendly, Menu Driven HMI Touch Screen Operator Interface
- Quick and Easy Tool-less Changeover
- Standard Laminar Airflow Design for Cleanroom
- Optional HEPA Filtration System Available
- Servo Motor Driven Technology

# Monobloc Specifications

## Additional Features:

- 💧 Handles free flowing to viscous products
- 💧 Fill Accuracy: +/- .5% or better
- 💧 Fill Volume Range: 0.5mL to 500mL
- 💧 Choice of Filling System: Positive Displacement, Peristaltic, or Micro-Dose Rotary Piston

## Standard MNB-2000 Specifications:



Filamatic offers **COMPLIMENTARY** product testing and purchasing options including rentals and lease arrangements. For more information call **866.258.1914** or email [info@filamatic.com](mailto:info@filamatic.com).

