

FILAMATIC®



**Specialists in Liquid Filling,
Capping and Labeling Solutions**

FlexFIL- Fully-Automated Solution

General Overview

- Unique ability to sort, feed, remove cap, fill and continue closing
- Choice of 2 -3 filling stations
- Optional HEPA filtration system available with laminar flow while still in a compact footprint
- User-friendly, menu driven HMI touchscreen interface
- Capable of running plastic and glass bottles, as well as different container shapes and sizes
- Servo-driven technology for improved control and reliability
- Quick and easy tool-less changeovers



Unique Performance. Undeniable Quality. Unrivalled Support.
Go Further, Do More with Filamatic.

***FlexFIL - Fully-Automated
Monobloc Filling Solution***

866.258.1914
info@filamatic.com
www.filamatic.com



Technical Features

- Remote diagnostics that are convenient and easy-to-access from across the globe
- Capable of bottom-up fill with locate fill and static fill profiles
- Choice of filling systems: Positive displacement, Peristaltic, or Micro-Dose Rotary Piston
- Ability to dispense a range of liquids: free-flowing, semi-viscous and heavy viscous products
- Fill accuracy: +/- 0.5% or better
- Optional gross or Net weigh capabilities
- 100% verification of cap/closing torque

PUMP SYSTEM	FILL STATIONS	VOLTS	CPM	FILL RANGE
Piston	2 - 3	208VAC	+ 25	Up to 1 liter
Peristaltic	2 - 3	208VAC	+ 25	Up to 1 liter

• CPM production will range with different container sizes, fill volumes, and varying product viscosities

FILAMATIC

4119 Fordleigh Road • Baltimore, MD 21215
866.258.1914 • info@filamatic.com • www.filamatic.com