



## DFS Series - Automatic Trolley Filling Systems

### General Overview

- Multiple frame size options available to handle a large number of pumps
- Dockable trolley system reduces changeover to +/- 5 minutes
- Tool-less adjustments
- Variety of metering systems available:
  - Flow meter pump
  - Lobe pump
  - Piston pump
- Bottom-up, locate, and static filling modes with nozzle position
- Fill rates of 10 to +200 CPM
- Fill accuracy +/- 0.5% or better
- Indexing and continuous motion filling



Unique Performance. Undeniable Quality. Unrivaled Support  
Go Further, Do More with FILAMATIC.



# DFS Series - Automatic Trolley Filling Systems

866.258.1914  
[info@filamatic.com](mailto:info@filamatic.com)  
[www.filamatic.com](http://www.filamatic.com)



## Technical Features

- Interlocked safety enclosure constructed of lexan and stainless steel
- Best for moderate to large fill volumes (50ml to Gallons)
- Comprehensive setup, usage, drawings, and validation documentation
- Fallen bottle, infeed bottle supply, discharge bottle back-up sensors on the container handling system
- Clean-in-place (CIP) or optional remote CIP system
- Recipe-driven system



FILAMATIC

4119 Fordleigh Road • Baltimore, MD 21215  
866.258.1914 • [info@filamatic.com](mailto:info@filamatic.com) • [www.filamatic.com](http://www.filamatic.com)