



## SmartFIL - Semi-Automatic Benchtops

### General Overview

- Easy-to-use HMI touchscreen
- Unmatched filling flexibility
- Self-priming and reversible for drip-free filling
- Recipe recall that allows you to create, save, and copy +100 recipes
- Stackable peristaltic pumps for increased production
- Advanced fill-weight compensation
- Cleanroom ready design
- Easy setup, low maintenance, and quick changeovers



Unique Performance. Undeniable Quality. Unrivaled Support.  
Go Further, Do More with FILAMATIC.


# SmartFIL - Benchtop Semi-Automatic Models

866.258.1914  
info@filamatic.com  
www.filamatic.com



### Technical Features

- Tool-less changeovers
- 110 V power
- Stainless steel cabinet
- Quick clean-up
- Fill accuracy +/- 0.5% depending upon metering pump selected
- Modular design with built-in expansion capability



Peristaltic	Fill Volume Range/Rates	Product Description
MasterFlex	.08 to 480ml Per Minute	Free Flowing Semi-Viscous
Watson Marlow	4.37 Liters Per Minute	Free Flowing Semi-Viscous
Flexicon	4.50 Liters Per Minute	Free Flowing Semi-Viscous

FILAMATIC

4119 Fordleigh Road • Baltimore, MD 21215  
 866.258.1914 • info@filamatic.com • www.filamatic.com